

# CLEANSTATION CS 750 – Clean, dry carts – every time

Resource-responsible, space-saving and process-efficient defines the Belimed CLEANSTATION CS 750 cart washer. It's the only real choice to meet today's demand for fast, effective and efficient cart and container processing in the CSSD.



*CLEANSTATION CS 750 washer, disinfectant and dryer for carts, containers, beds and other large-volume washing items.*

## **Resource-responsible**

A highly efficient tank wash system, coupled with an on-demand rinse water heating design delivers lowest-in-class water use. This unique combination also eliminates the need to cool drain discharge water, saving even more of the valuable natural resource, water, as well as the financial resource-operational cost.

## **Space-saving**

The tinted glass pivot doors, compact wash tank and plumbing design all result in a reduced footprint while increasing the chamber volume, returning valuable floor space for use for additional critical needs of the department.

## **Process-efficient**

Our high-output dryer, tilting floor and high-temperature rinse all ensure that the most time-consuming step, drying, is reduced to an absolute minimum. Clean and dry carts in 6 to 7 minutes, with no compromise – unprecedented throughput so the packing area isn't left waiting.

The large, bright status indicator, easily visible from across the room, lets personnel know that the cycle is complete and ready to be unloaded, further improving productivity.

**Belimed**  
Infection Control

System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed, Inc., 2284 Clements Ferry Road, Charleston, SC 29492, Phone 843 216 7424, Fax 843 216 7707, info@belimed.us, [www.belimed.us](http://www.belimed.us)



Status of the machine can be quickly assessed through the large, bright indicator, easily seen from across the room.



The intuitive user interface, with illuminated touch-sensitive buttons, ensures proper operation.

### User-friendly

A fully walkable floor makes cart loading and retrieval easier while the intuitive color display provides the user with a clear, simple interface to all functions. With sleek glass doors on both sides and internal illumination, for full visibility into the chamber, the CS 750 was designed with the user in mind.

### Flexible design

Available with three chamber lengths, short, medium or long, as well as with a left- or right-hand service area, the CS 750 can be configured to fit within your available space while also meeting your processing needs.

## Technical Data

Parameter	CS 750S	CS 750M	CS 750L
Usable chamber dimensions H x W x D in (mm)	78.7" x 39.4" x 53.1" (2000 x 1000 x 1350)	78.7" x 39.4" x 90.6" (2000 x 1000 x 2300)	78.7" x 39.4" x 116.1" (2000 x 1000 x 2950)
Outer dimensions H x W x D in (mm)	118.5" x 98.4" x 88.6" (3010 x 2500 x 2250)	118.5" x 98.4" x 126.0" (3010 x 2500 x 3200)	118.5" x 98.4" x 151.6" (3010 x 2500 x 3850)
Pit depth in (mm)	4.7" (120)	4.7" (120)	4.7" (120)



System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed, Inc., 2284 Clements Ferry Road, Charleston, SC 29492, Phone 843 216 7424, Fax 843 216 7707, info@belimed.us, www.belimed.us