



Water Treatment for Sterile Processing



- Compact Footprint
- Simple, Integrated System
- Easy To Operate
- Exceptional Performance
- Competitive Price



Water Treatment Solutions for Sterile Processing

Complete water treatment systems for central sterile services from AmeriWater are designed to use a minimum of valuable floor space and provide high quality water for washing, rinsing and disinfecting. Field-proven in medical facilities across North America, AmeriWater's Sterile Processing water treatment systems are designed to minimize maintenance and deliver greater energy efficiency while providing cleaner water for cleaning, disinfecting, and sterilization of medical equipment.

Combination RO/DI for the Sterile Processing Department

This complete system takes up little floor space, yet includes everything you need to feed purified water to all of your washers, disinfectors and sterilizers in the CSS department. This system includes:

- Prefiltration for sediment removal
- Onboard anti-scalant to reduce hardness
- Reverse Osmosis to purify your water
- Deionizer to polish the RO water
- Storage tank with distribution pump to assure continuous pure water to your systems

Local installation and service available.



Exchangeable Resin Pack Deionizer for Single Applications

The Silex system is a dual PVC column designed with a 200,000 ohm water quality monitor to indicate the system is producing purified water for a single washer/disinfecter or sterilizer with very little maintenance. Recirculation pump keeps the water moving, which helps maintain high quality water at all times. When the indicator light changes from green to red (and alarm sounds), simply remove the lids and exchange the resin packs. Resin packs are not disposed of but regenerated by AmeriWater, which helps keep your CSS department green.



CERTIFICATIONS

ISO 13485/FDA 1530185
ETL ANSI/UL 61010
Seismic OPA-2202-07

Custom configured as part of the complete Belimed system solution

Belimed
Infection Control

Belimed, Inc.
2284 Clements Ferry Road
Charleston, SC 29492
Tel: 800-457-4117
Fax: 843-216-7707
www.belimed.us