

# BLANKET & FLUID WARMING CABINETS

MAC Medical, Inc. offers a variety of Blanket and Fluid Warming Cabinets. As with all of our products, warming cabinets undergo stringent testing procedures during the manufacturing process to assure a long life, the highest quality and accurate temperature settings. We offer standard sizes along with custom sizes to meet specific requirements.

## Standard Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2 - 6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and visual over-temp alarm
- UL/CUL approved
- FDA registered
- OSHPD pre-approved

## Options

- Full view glass door
- Mobile dolly or cart
- 18" depth
- Full recess
- Pass through chamber
- Left hand door hinge
- Keyed door lock
- Temperature recorder
- Seismic brackets
- 220/240 V, 50/60HZ
- Custom sizes



(Model TWC24-G-TL-MB shown above on optional dolly)

### Model # TWC24

**Overall Size**  
30"W x 26.5"D x 78.25"T  
14.00 cubic feet of storage  
1 adjustable shelf

**Upper Chamber I.D.**  
26"W x 23"D x 11.25"T  
3.9 cubic feet of storage  
Temperature range: 90-110 F.  
Accuracy +/- 1 degree F. of set point

**Intermediate Chamber**  
26"W x 23"D x 10.75"T  
3.7 cubic feet of storage  
Temperature range: 90-135 F.  
Accuracy +/- 1 degree F. of set point

**Lower Chamber I.D.**  
26"W x 23"D x 18.50"T  
6.4 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

**Electrical**  
120V, 13.5 amp., 50/60HZ  
.65 kw/hr, 2238 btu/hr  
UL/CUL listed  
FDA registered



(Model DWC24-TL-MB shown above on optional dolly)

### Model # DWC24

**Overall Size**  
30"W x 26.5"D x 74.5"T  
17.2 cubic feet of storage  
3 adjustable shelves

**Upper Chamber I.D.**  
26"W x 23"D x 15.25"T  
5.27 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

**Lower Chamber I.D.**  
26"W x 23"D x 34.5"T  
11.93 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 3 degree F. of set point

**Electrical**  
120V, 12 amp., 50/60HZ  
.75 kw/hr, 2557 btu/hr  
UL/CUL listed  
FDA registered



(Model DWC2424-TL-MB shown above on optional dolly)

### Model # DWC2424

**Overall Size**

24"W x 26.5"D x 74.5"T  
15.95 cubic feet of storage  
3 adjustable shelves

**Upper Chamber I.D.**

20"W x 23"D x 15.25"T  
4.05 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

**Lower Chamber I.D.**

26"W x 23"D x 34.5"T  
11.9 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 3 degree F. of set point

**Electrical**

120V, 10 amp., 50/60HZ  
.46 kw/hr, 1570 btu/hr  
UL/CUL listed  
FDA registered



(Model DWC2460-G-TL shown above)

### Model # DWC2460

**Overall Size**

24"W x 20.5"D x 64.75"T  
7.45 cubic feet of storage  
3 adjustable shelves

**Upper Chamber I.D.**

20"W x 17"D x 13.5"T  
2.65 cubic feet of storage  
Temperature range: 90-160 F  
Accuracy +/- 1 degree F. of set point

**Lower Chamber I.D.**

20"W x 17"D x 24.5"T  
4.80 cubic feet of storage  
Temperature range: 90-160 F  
Accuracy +/- 2 degree F. of set point

**Electrical**

120V, 10 amp., 50/60HZ  
.49 kw/hr, 1678 btu/hr  
UL/CUL listed  
FDA registered



(Model DWC36-G-TL-MB shown above on optional dolly)

### Model # DWC36

**Overall Size**

30"W x 26.5"D x 74.5"T  
17.59 cubic feet of storage  
4 adjustable shelves

**Upper Chamber I.D.**

26"W x 23"D x 25"T  
8.63 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 2 degree F. of set point

**Lower Chamber I.D.**

26"W x 23"D x 26"T  
8.96 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 2 degree F. of set point

**Electrical**

120V, 12 amp., 50/60HZ  
.69 kw/hr, 2356 btu/hr  
UL/CUL listed  
FDA registered



(Model DWC60-G-TL shown above)

### Model # DWC60

**Overall Size**

30"W x 20.5"D x 64.75"T  
9.73 cubic feet of storage  
3 adjustable shelves

**Upper Chamber I.D.**

26"W x 17"D x 13.5"T  
3.46 cubic feet of storage  
Temperature range: 90-160 F  
Accuracy +/- 1 degree F. of set point

**Lower Chamber I.D.**

26"W x 17"D x 24.5"T  
6.27 cubic feet of storage  
Temperature range: 90-160 F  
Accuracy +/- 2 degree F. of set point

**Electrical**

120V, 10 amp., 50/60HZ  
.49 kw/hr, 1678 btu/hr  
UL/CUL listed  
FDA registered



**Model # SWC1518**

**Overall Size**  
18"W x 15.5"D x 22"T

**Chamber I.D.**  
14"W x 13"D x 13"T  
1.4 cubic feet of storage  
Temperature range: 90-135 F.  
Accuracy +/- 1 degree F. of set point

**Electrical**  
120V, 3.5 amp., 50/60HZ  
.39 kw/hr, 1332 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC1518-G-TL shown above)



**Model # SWC1824**

**Overall Size**  
24"W x 20.5"D x 24.5"T

**Chamber I.D.**  
20"W x 17"D x 15.25"T  
3.0 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

**Electrical**  
120V, 3.5 amp., 50/60HZ  
.45 kw/hr, 1537 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC1824-G-TL shown above)



**Model # SWC24**

**Overall Size**  
30"W x 26.5"D x 24.5"T  
1 adjustable shelf

**Chamber I.D.**  
26"W x 23"D x 15.25"T  
5.27 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 1 degree F. of set point

**Electrical**  
120V, 7.5 amp., 50/60HZ  
.30 kw/hr, 1024 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC24-G-TL shown above)



**Model # SWC2460**

**Overall Size**  
24"W x 20.5"D x 64.75"T  
3 adjustable shelves

**Chamber I.D.**  
20"W x 17"D x 49.5"T  
9.73 cubic feet of storage  
Temperature range: 90-160 F  
Accuracy +/- 1 degree F. of set point

**Electrical**  
120V, 7.5 amp., 50/60HZ  
.49 kw/hr, 1678 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC2460-G-TL-MB shown above on optional dolly)



**Model # SWC36**

**Overall Size**  
30"W x 26.5"D x 36"T  
2 adjustable shelves

**Chamber I.D.**  
26"W x 23"D x 25"T  
8.63 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 2 F. degree of set point

**Electrical**  
120V, 7.5 amp., 50/60HZ  
.35 kw/hr, 1195 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC36-G-TL shown above)



**Model # SWC60**

**Overall Size**  
30"W x 20.5"D x 64.5"T  
3 adjustable shelves

**Chamber I.D.**  
26"W x 17"D x 51"T  
13.05 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 3 degree F. of set point

**Electrical**  
120V, 7.5 amp., 50/60HZ  
.46 kw/hr, 1570 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC60-G-TL-MB shown above on optional dolly)



**Model # SWC72**

**Overall Size**  
30"W x 26.5"D x 74.5"T  
2 adjustable shelves  
2 fixed shelves

**Chamber I.D.**  
26"W x 23"D x 61"T  
21 cubic feet of storage  
Temperature range: 90-160 F.  
Accuracy +/- 3 degree F. of set point

**Electrical**  
120V, 7.5 amp., 50/60HZ  
.53 kw/hr, 1810 btu/hr  
UL/CUL listed  
FDA registered

(Model SWC72-G-TL-MB shown above on optional dolly)

**Optional Mobile Stands:**



Half size and counter top units  
(If ordered with warmer add -MS to the model #)

**T0007:** 26"W x 34"D x 34"T  
Recommended for models SWC24 and SWC36

**T0023:** 26"W x 20"D x 34"T  
Recommended for models SWC1518 and SWC1824

**Optional Mobile Base:**



Newly Redesigned Mobile Base:

**T0087:** 33"W x 28"D x 6.75"T  
(If ordered with the warmer, add -MB to the model #.  
Note this will add 2" to the standard warmer height)



Standard keyed temperature lockout



Optional door lock

## How to order Blanket & Fluid Warming Cabinets

Options Available	18" Depth	Glass Doors	Mobile Base	Mobile Stand	Recessed Unit	LH Door Swing, Hinges on Left	Keyed Door Lock	Temperature Recorder	220/240 V, 50/60 HZ	Pass Through Chamber
Option Extension	18	G	MB	MS	R	LH	DL	TR	220	P

For Example: Model # DWC24 with the 18" depth option, Glass Doors, Recessed unit, Hinges on the left, Keyed door lock, 220/240 V 50/60 HZ = Order Model # DWC24-18-G-R-LH-DL-220

If you require any assistance please call us with any questions at (618) 476-3550 or toll free at 1-877-828-9975.

## Standard Models

Model #	Outside Dimensions
<b>Triple Chamber Units</b>	
TWC24	30"W x 26.5"D x 78.25"T
TWC24-G	30"W x 26.5"D x 78.25"T
<b>Double Chamber Units</b>	
DWC24	30"W x 26.5"D x 74.5"T
DWC24-G	30"W x 26.5"D x 74.5"T
DWC2424	24"W x 26.5"D x 74.5"T
DWC2424-G	24"W x 26.5"D x 74.5"T
DWC2460	24"W x 20.5"D x 64.75"T
DWC2460-G	24"W x 20.5"D x 64.75"T
DWC36	30"W x 26.5"D x 74.5"T
DWC36-G	30"W x 26.5"D x 74.5"T
DWC60	30"W x 20.5"D x 64.75"T
DWC60-G	30"W x 20.5"D x 64.75"T
<b>Single Chamber Units</b>	
SWC1518	18"W x 15.5"D x 22"T
SWC1518-G	18"W x 15.5"D x 22"T
SWC1824	24"W x 20.5"D x 24.5"T
SWC1824-G	24"W x 20.5"D x 24.5"T
SWC24	30"W x 26.5"D x 24.5"T
SWC24-G	30"W x 26.5"D x 24.5"T
SWC2460	24"W x 20.5"D x 64.75"T
SWC2460-G	24"W x 20.5"D x 64.75"T
SWC36	30"W x 26.5"D x 36"T
SWC36-G	30"W x 26.5"D x 36"T
SWC60	30"W x 20.5"D x 64.75"T
SWC60-G	30"W x 20.5"D x 64.75"T
SWC72	30"W x 26.5"D x 74.5"T
SWC72-G	30"W x 26.5"D x 74.5"T

**NOTE:** Due to the vast number of sizes we have available, not all of our models are shown. Please let us know your specific dimensions to get the model that will best suit your needs. **If we don't currently have the size you need, we will build a custom one for you!**

## WARRANTY

All warming cabinets manufactured by MAC Medical, Inc. will carry a lifetime guarantee against product craftsmanship, 1-year labor and 2-year parts guarantee. The factory will service all units without cost to the buyer for 1-year from shipment. After the 1-year period, replacement of a defective part (labor) will be at buyer's expense. We will exchange all defective parts at no cost to the buyer for a period of 2-years from shipment. All defective parts must be returned within 30 days to ensure proper credit. An RMA from MAC Medical, Inc. must be obtained prior to items return.