

BELICHEM 100

BELICHEM 100 is a liquid alkaline detergent for use in ultrasonic cleaners and automated washer machines for the removal of heavy contamination on surgical instruments, anesthetic utensils, and in laboratories.

Contains potassium silicate, potassium hydroxide & chelants and is compatible with stainless steel. This product is not compatible with anodized aluminum and most non-ferrous metals.

FEATURES

Contains chelating agents
Contains surfactants
Highly concentrated
Low foaming
Biodegradable

BENEFITS

Hard water tolerant
Optimal wetting and detergency
Low cost per use
Effective with all automated washers
Non-polluting, ecology minded

DIRECTIONS FOR USE

Use 3.5 ml to 5 ml per liter of water or 1/3 ounce to 2/3 ounce per gallon of water. Dilutions may vary due to water conditions, soil load, and type of equipment. Contact your Belimed Sales Consultant for recommended usages.

PHYSICAL PROPERTIES

Appearance:	Clear, straw colored liquid
Fragrance:	None
pH 1% @ 20°C:	11.6
SG @ 20°C	1.17
Foam:	Negligible
Available Alkalinity:	4.7

Available in 5, 15, 30, and 55 gallon containers